


# AIRWAY MANAGEMENT

---





# MUCUS SPECIMEN TRAPS






MST-3070



MST-3000



MST-4000



Mucus specimen traps designed to obtain mucus specimen for microbiological examination.

- Sterile, individual pack
- 25ml or 70ml trap
- Attachment for suction catheter and suction connecting tube
- Patient details label
- Spare cap, screw fit

Product Code	Capacity	Type	Qty Per Box
MST-3000	25ml	With spare screw cap - 14cm tubing	50
MST-3070	25ml	With reduced length connection to suction catheter to ease movement of secretions to the trap. With spare screw cap & patient label.	50
MST-3180	25ml	With spare screw cap - 29cm tubing	75
MST-4000	70ml	For use in bronchoscopy. - 11cm tubing. With spare screw cap	50
MST-4000/S	70ml	For use in bronchoscopy. - 1mt tubing. With spare screw cap	25
MST-5000	25ml	With spare screw cap - 9cm tubing with VC	50
MST-7000	70ml	With male connector	50



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)



---

Ivor Shaw Limited t/a Pennine Healthcare  
300 City Gate, London Road, Derby, DE24 8WY, UK.  
Tel: +44 (0)1332 794880  
[www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)

 @arma\_health  Pennine Healthcare  @Pennine\_Health

BRO-5\_Rev\_4 02/2023