

Pressure Infusion Cuffs



REF: 43.018.00.010, 43.019.00.010

Product group: Infusion/Transfusion Accessories

Intended use:

The pressure infusion cuffs are designed for the rapid infusion of fluids and for use with various invasive pressure monitoring systems.

Product description:

For rapid fluid infusion and use with invasive pressure monitoring systems. Features an inflatable cuff with transparent front, pressure gauge, three-way stopcock, and inflation bulb.

Not MRI compatible.

Sterile: No

REF	Capacity	Selling unit quantity	Selling unit GTIN	Packing unit quantity	Packing unit GTIN
43.018.00.010	500 ml	1 pcs	04034342375805	50 pcs	04034342376314
43.019.00.010	1000 ml	1 pcs	04034342375812	50 pcs	04034342376321

Materials

Product	Nylon, PVC, ABS, PC, POM, steel
Packaging	PE bag, cardboard, corrugated cardboard
Latex	No
DEHP/BBP/DBP	No

Additional information

Regulatory basis	Regulation (EU) 2017/745
Notified Body	-
Classification	I
UMDNS	13-100
eClass	34221219
Manufacturer	P. J. Dahlhausen & Co. GmbH
Single-use Product	Yes
Sterilisation	-
Shelf life	2 Years
Instructions for Use	Yes