

## Silicone balloon catheter, Foley type


**REF:**

23.000.14.208, 23.000.14.210, 23.000.14.212,  
23.000.14.214, 23.000.14.216, 23.000.14.218,  
23.000.14.220, 23.000.14.222, 23.000.14.224,  
23.000.24.218, 23.000.24.220, 23.000.24.222,  
23.000.24.224

**Product group:**

Balloon Catheters

**Intended use:**

The balloon catheter is used to drain urine from the bladder and then transfers it to a urine collection container.

**Product description:**

The silicone balloon catheter type Foley consists of a 2-lumen shaft with a proximal funnel and a valve for the distal retention balloon.

REF: 23.100.14.208 and 23.100.14.210 have a guide wire as an insertion aid and are 30 cm long.

From REF 23.000.14.212 the length is 40 cm.

<b>Allergenic potential:</b>	Latex-free
<b>Eye arrangement:</b>	Oposing
<b>Application area:</b>	Male/female
<b>Number (eyes):</b>	2
<b>Number (lumen):</b>	2-lumen
<b>X-ray contrastable:</b>	X-ray contrast strip
<b>Sterile:</b>	Yes

REF	Outer diameter	Balloon capacity	Colour-coding	Selling unit quantity	Selling unit GTIN	Packing unit quantity	Packing unit GTIN
23.000.14.208	8 Ch	3 ml	Black	10 pcs	04034342081843	500 pcs	04034342328382
23.000.14.210	10 Ch	3 ml	Grey	10 pcs	04034342081874	500 pcs	04034342328399
23.000.14.212	12 Ch	5–10 ml	White	10 pcs	04034342081904	500 pcs	04034342328405
23.000.14.214	14 Ch	5–10 ml	Green	10 pcs	04034342081935	500 pcs	04034342328412
23.000.14.216	16 Ch	5–10 ml	Orange	10 pcs	04034342081966	500 pcs	04034342328429
23.000.14.218	18 Ch	5–10 ml	Red	10 pcs	04034342081997	500 pcs	04034342328436
23.000.14.220	20 Ch	5–10 ml	Yellow	10 pcs	04034342082024	500 pcs	04034342328443
23.000.14.222	22 Ch	5–10 ml	Violet	10 pcs	04034342082055	500 pcs	04034342328450
23.000.14.224	24 Ch	5–10 ml	Blue	10 pcs	04034342082086	500 pcs	04034342328467
23.000.24.218	18 Ch	20 ml	Red	10 pcs	04034342078355	500 pcs	04034342328481
23.000.24.220	20 Ch	30 ml	Yellow	10 pcs	04034342078386	500 pcs	04034342328498
23.000.24.222	22 Ch	30 ml	Violet	10 pcs	04034342082215	500 pcs	04034342328504
23.000.24.224	24 Ch	30 ml	Blue	10 pcs	04034342082246	500 pcs	04034342328511

## Materials

Product	Silicone, ABS
Packaging	PP/PE, paper, cardboard, corrugated cardboard
Latex	No
DEHP/BBP/DBP	No

## Additional information

Regulatory basis	Directive 93/42/EEC
Notified Body	TÜV SÜD Product Service GmbH (0123)
Classification	IIb
UMDNS	10-720
eClass	34250301
Manufacturer	P. J. Dahlhausen & Co. GmbH
Single-use Product	Yes
Sterilisation	EO gas
Shelf life	5 Years
Instructions for Use	Yes